

EEG Acquisition (ACK) Block

- (1) Start SSVEP block in background with $(f1, f2)$ as frequency parameters
- (2) Start SP block in background
- (3) Create a socket connection between this block and SP block
- (4) Initialize Emotiv EEG
- (5) Execute main BCI loop for n_trials times (or until interrupted)
 - (a) Acquire resting EEG for *duration* seconds
 - (b) Send signal to SSVEP block to start flickering
 - (c) Acquire SSVEP EEG for *duration* seconds
 - (d) Send signal to SSVEP block to stop flickering
 - (e) Send EEG data to SP block for processing over the socket

