EEG Acquisition (ACK) Block

(1) Start SSVEP block in background with (f_1, f_2) as frequency parameters (2) Start SP block in background (3) Create a socket connection between this block and SP block (4) Initialize Emotiv EEG (5) Execute main BCI loop for *n_trials* times (or infinitely) (a) Acquire resting EEG for duration seconds (b) Send signal to SSVEP block to start flickering (c) Acquire SSVEP EEG for duration seconds (d) Send signal to SSVEP block to stop flickering (e) Send EEG data to SP block for processing over the socket Socket Block until the arrival of a signal Listen the socket for new EEG dataSignal received 🖣 Start flickering the LEDs Process EEG data and select between **left** (f_1) or **right** (f_2).Signal received Signal Processing (SP) Block SSVEP Stimulation Block